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Alarms and Alerts

Table 7-1: Alarm and Alert Conditions (continued)

Function	Alert icon	Response	Do this
Single channel leak	Leak indicator icon	The ventillator has detected a leak. Indicators will be flashing and the ventillator will not breathe.	1. Alerted to the patient first, as instructed by your clinician. 2. Check the patient circuit and accessory equipment. Make sure that all equipment is seated, connected and properly secured. Replace if necessary. 3. If the ventillator will not breathe, disconnect the patient circuit and check for a leak in the patient circuit and ventillator tubing. 4. If the leak persists, disconnect the patient circuit and contact your clinician for guidance.
Single alarm tone	Alarm icon	The ventillator has detected a non-critical alarm.	1. Alerted to the patient first, as instructed by your clinician. 2. Contact your clinician. Do not use the ventillator until a warning message is received.
Residual oxygen lamp	Residual oxygen lamp icon	The ventillator has detected a low residual oxygen alarm.	1. Acknowledge the alarm indicator by pressing Alarm Silence/Reset. 2. If the ventillator has switched to a normal battery, ensure ventilation is provided as normal power source as soon as possible.
None	None	The ventillator has no alarm or error.	1. Alerted to the patient first, as instructed by your clinician. 2. Check the ventilation settings. Make sure the ventilation settings are set according to the prescription. 3. If the alarm condition persists, disconnect use of the ventillator and contact your clinician for guidance.

- 7.2 Resetting Alarms
- After you have corrected the alarm condition, press Alarm Silence/Reset to reset the alarms and deactivate the indicators on the Alarm panel.
- 7.3 Alarm Latching
- Some of the ventillator's low pressure alarms can operate in either a latching or non-latching mode.
- A non-latching alarm does not require you to press Alarm Silence/Reset to stop the audible alarm when the condition has been corrected. It will automatically reset if the audible alarm has the visual indicator on the Alarm panel will still light. This lets you know that an alarm event occurred and has corrected.
 - A latching alarm will not stop the audible alarm unless the condition has been corrected and the Alarm Silence/Reset key has been pressed. Pressing the Alarm Silence/Reset key on the ventillator does not stop the latching mode of the alarm.
- To stop the latching mode of the alarm:
- Press Mute/EC.
 - Press the UP Arrow or Down Arrow keys until the display screen reads: Press START/ENTER to CHANGE ALARM LATCHING STATUS.

Alarms and Alerts

Table 7-1: Alarm and Alert Conditions (continued)

Function	Alert icon	Response	Do this
Pre-Silent Audible Alarms	Pre-silent audible alarm icon	There may be situations in which you wish to silence an anticipated alarm for a brief period. The ventillator allows you to "pre-silent" audible alarms for up to eight seconds.	1. To pre-silent audible alarms, press Alarm Silence/Reset. This will silence any audible alarms for up to eight seconds. During this period, the indicators on the alarm display panel will not light to indicate an alarm condition, but the alarm will not sound. 2. To stop pre-silencing of audible alarms, press Alarm Silence/Reset.

- 7.4 Pre-Silencing Audible Alarms
- There may be situations in which you wish to silence an anticipated alarm for a brief period. The ventillator allows you to "pre-silent" audible alarms for up to eight seconds.
- To pre-silent audible alarms, press Alarm Silence/Reset. This will silence any audible alarms for up to eight seconds. During this period, the indicators on the alarm display panel will not light to indicate an alarm condition, but the alarm will not sound.
 - To stop pre-silencing of audible alarms, press Alarm Silence/Reset.
- If an alarm condition occurs and is corrected during the pre-silence period, pressing Alarm Silence/Reset will reset the alarm and deactivate the indicators.

Alarms and Alerts

Table 7-1: Alarm and Alert Conditions (continued)

Function	Alert icon	Response	Do this
Regulating level of alarm panels	Regulating level of alarm panels icon	The pressure in the patient circuit has dropped below the preset alarm level.	1. Alerted to the patient first, as instructed by your clinician. The patient circuit is in danger: continue to breathe. 2. Check the patient circuit for leaks or loose connections. Fix or replace the patient circuit if necessary. 3. Inspect for and remove water from small tubing. 4. Check the ventillator settings. Verify that the ventillator settings are set according to the prescription. 5. If the alarm condition persists, disconnect use of the ventillator and contact your clinician for guidance.
Alarm Silence/Reset indicator	Alarm Silence/Reset indicator icon	The ventillator has detected an alarm.	1. Alerted to the patient first, as instructed by your clinician. 2. Check the patient circuit for leaks or loose connections. Fix or replace the patient circuit if necessary. 3. Press the Silence key (S). If this corrects the alarm, ensure normal ventilation. If the alarm has not cleared, press the Stop key (O). 4. Using the ventillator from AC power and read 30 seconds. Plug the ventillator into AC power and ensure ventilation. If the alarm has not cleared, proceed to step 3. 5. Using the ventillator from AC power, this panel will flash (DBU) for 30 seconds or more seconds. Plug the ventillator into AC power and ensure ventilation. If the alarm has not cleared, contact your clinician.
Normal battery charge is depleted	Normal battery charge is depleted icon	The normal battery charge is depleted.	1. Disconnect the ventillator from the patient. 2. Press and hold the Silence key (S) for three (3) seconds. 3. Contact the ventillator for an external power supply. If using AC power, the normal battery will automatically begin to charge. 4. Press Mute/EC to resume ventilation. 5. If using a charged normal battery, you must press Silence/Reset and reset (S) your ventillator. The battery will continue to charge. 6. Operate the ventillator on AC power for at least four (4) hours to recharge the normal battery.

Alarms and Alerts

Table 7-1: Alarm and Alert Conditions (continued)

Function	Alert icon	Response	Do this
High flow alarm	High flow alarm icon	The ventillator has detected a high flow alarm.	1. Alerted to the patient first, as instructed by your clinician. 2. Check the patient circuit for leaks or loose connections. Fix or replace the patient circuit if necessary. 3. Check the ventillator settings and verify they are set according to the clinician's prescription. 4. If the alarm condition persists, disconnect use of the ventillator and contact your clinician for guidance.
High flow alarm (PF)	High flow alarm (PF) icon	The ventillator has detected a high flow alarm (PF).	1. Alerted to the patient first, as instructed by your clinician. 2. Check the patient circuit for leaks or loose connections. Fix or replace the patient circuit if necessary. 3. Check the ventillator settings and verify they are set according to the clinician's prescription. 4. If the alarm condition persists, disconnect use of the ventillator and contact your clinician for guidance.
High flow alarm (PF) (DBU)	High flow alarm (PF) (DBU) icon	The ventillator has detected a high flow alarm (PF) (DBU).	1. Alerted to the patient first, as instructed by your clinician. 2. Check the patient circuit for leaks or loose connections. Fix or replace the patient circuit if necessary. 3. Check the ventillator settings and verify they are set according to the clinician's prescription. 4. If the alarm condition persists, disconnect use of the ventillator and contact your clinician for guidance.
Single beep, repeat every 5 minutes	Single beep, repeat every 5 minutes icon	There is a condition with use or status of the wiring harness.	1. Alerted to the patient first, as instructed by your clinician. 2. Check the ventillator settings and verify they are set according to the clinician's prescription. 3. If the alarm condition persists, disconnect use of the ventillator and contact your clinician for guidance.
Single beep, repeat every 30 minutes	Single beep, repeat every 30 minutes icon	The normal battery charge is depleted.	1. Disconnect the ventillator from the patient. 2. Press and hold the Silence key (S) for three (3) seconds. 3. Contact the ventillator for an external power supply. If using AC power, the normal battery will automatically begin to charge. 4. Press Mute/EC to resume ventilation. 5. If using a charged normal battery, you must press Silence/Reset and reset (S) your ventillator. The battery will continue to charge. 6. Operate the ventillator on AC power for at least four (4) hours to recharge the normal battery.
Single beep, repeat every 15 minutes	Single beep, repeat every 15 minutes icon	The ventillator has detected a low battery level.	1. Alerted to the patient first, as instructed by your clinician. 2. Check the patient circuit and accessory equipment. Make sure that all equipment is seated, connected and properly secured. Replace if necessary. 3. If the ventillator will not start ventilation, disconnect use of the ventillator and contact your clinician for guidance.

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